

Medical mouth & nose protection



Filtration half mask - respirators



Cloth face covering – miamasuku®



Purpose	Protection of strangers	Self-protection	Protection of strangers
Fit	– loose-fitting – prevent large-particle droplets	– tight fitting – reduces the wearer’s exposure to respiratory contaminants	– loose-fitting – reduces the spread of droplets
Period of use	Disposable product Dispose of after single use. At the latest as soon as they have become damp.	Disposable product Change after max 8 hours of use.	Zero waste Unlimited reusable!
Classification	Medical device with CE certificate and DIN EN 14683:2019-6	Medical device with CE certificate and standard DIN EN 149:2001-10	Hygiene product